

Biological and Chemical Indicator

Sterility Assurance Test Records for monitoring the System 1 sterilization process

| Date | Time | Processor Serial Number and/or Cycle # | Location and/or Load | Operator | Incubator Temp | Lot Number SS (spore strips) CM(culture media) | Results | | | CI ✓ ✗ | Notes | Date/Time Out Initials |
|------|------|--|----------------------------|----------|-------------------|--|---------|----------|-------|-----------|-------|---------------------------|
| | | | | | | | Control | Test | | | | |
| | | | | | | | | 24 hours | Final | | | |
| | | | | | | S.S. | + - | + - | + - | | | |
| | | | | | | C.M. | | | | | | |
| | | | | | | S.S. | + - | + - | + - | | | |
| | | | | | | C.M. | | | | | | |
| | | | | | | S.S. | + - | + - | + - | | | |
| | | | | | | C.M. | | | | | | |
| | | | | | | S.S. | + - | + - | + - | | | |
| | | | | | | C.M. | | | | | | |
| | | | | | | S.S. | + - | + - | + - | | | |
| | | | | | | C.M. | | | | | | |
| | | | | | | S.S. | + - | + - | + - | | | |
| | | | | | | C.M. | | | | | | |
| | | | | | | S.S. | + - | + - | + - | | | |
| | | | | | | C.M. | | | | | | |
| | | | | | | S.S. | + - | + - | + - | | | |
| | | | | | | C.M. | | | | | | |
| | | | | | | S.S. | + - | + - | + - | | | |
| | | | | | | C.M. | | | | | | |
| | | | | | | S.S. | + - | + - | + - | | | |
| | | | | | | C.M. | | | | | | |

Biological

- + Turbidity and/or color change to or towards yellow
- No turbidity and no color change (sterilization conditions achieved)

Chemical

- ✓ Color Loss Complete (sterilant potency correct)
- ✗ Color Loss Incomplete

